



FACILITIES  
Emtech Laboratories  
7745 Garland Circle  
Roanoke, VA 24019

MAILING  
Emtech Laboratories  
P.O. Box 12900  
Roanoke, VA 24022

PHONE  
Local: (540) 265-9156  
Toll Free: (800) 336-5719  
Fax: (540) 265-9164

## Anesthesiologist & Stethoscope Earpieces

An anesthesiologist or stethoscope mold is a receiver with a ring. Attached is a metal adaptascope with 1 ft. of tubing, more if requested.

We recommend a vent for maximum comfort. Either can be made in any material. Lucite is the most popular. If requested we can make a stethoscope mold a full shell with an open sound channel for the stethoscope tip.

We can also make this mold with a sound bore and tubing for use with a hearing aid.

