



**FACILITIES**  
Emtech Laboratories  
7745 Garland Circle  
Roanoke, VA 24019

**MAILING**  
Emtech Laboratories  
P.O. Box 12900  
Roanoke, VA 24022

**PHONE**  
Local: (540) 265-9156  
Toll Free: (800) 336-5719  
Fax: (540) 265-9164

## The Janssen

The Janssen earpiece leaves the canal and entrance unobstructed. The deeper the sound outlet, the greater the usable gain. An increase of 8 to 10 db in gain before feedback can be obtained as the sound outlet is moved from the entrance to the canal inward 14 to 16 mm. The canal length and bore should be specified. This mold is available in Lucite.

